

Shoprider -STREAMER

ASK FOR THE BRAND SHOPRIDER



STREAMER

SHOPRIDER®
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

STREAMER

SHOPRIDER
www.shoprider.com.au
mobility

The SHOPRIDER™ **Streamer** is another great powerchair from Shoprider. This budget model offers many of the features found on much more expensive powerchairs. The Streamer comes with a high back adjustable seat, adjustable width & height armrests, flat free tyres and finger control steering. The Streamer is rear wheel drive and combined with two powerful 0.8hp motors it will deliver outstanding performance. This powerchair offers excellent value for money and reliability only a Shoprider can offer.

| SPECIFICATIONS | Streamer | FEATURES: |
|--------------------------------|------------|--------------------------------------|
| Length | 966mm | High back seat |
| Width | 560mm | Sleek body design |
| Height | 1016mm | Smooth operation |
| Battery Capacity (Minimum Amp) | 2 x 26Amp | Variable seat height |
| #Maximum Range | up to 13km | Adjustable arm rests |
| Maximum Load Single User* | 136Kgs | Controller interchangeable both arms |
| Charger off board | 5Amp | Free wheel device |
| Motor Size | 0.8 Hp | Option of swing away leg rest |
| Maximum Speed | 7.5Kph | Option of elevating leg rest |
| *Safe Climbing Angle | 9 Degrees | Option of Stump Support |
| *Max Climbing Angle | 9 Degrees | |
| Turning Radius | 650mm | |
| Wheel Size: | | |
| - Front Caster | 2 x 150mm | |
| - Rear drive wheel | 2 x 265mm | |

Driver weight may exceed weight of unit. Reduce speed when turning

IMPORTANT

Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



Streamer

9' (15%) -91kg
6' (10%) -113kg
0' (0) -136kg

SHOPRIDER
www.shoprider.com.au
mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au