

REDGUM BRAND
Daily Living and
Lifestyle Products

COMFORT LUX

BEDS AND MATTRESSES RANGE



*Compliant to:

- AS/NZS IEC 60601.1:2015
- IEC 60601-2-52 (Edition 1.1 2015-03)


COMFORT
QUALITY
VALUE

ASK FOR THE BRAND

redgumbrand.com.au

Comfort LUX - RG66BHL

High - Low King Size Single Bed



Upholstery

FUNCTIONS



ANTI-SNORE:

The adjustable bed head segment rises to an approximate 15 degree angle. This small movement can open up partially restricted airways.

FEATURES

- Multifunction - "Zero - G" and "Antisnore"
- Battery backup
- High - Low platform height that improves user safety
- Locking castors included
- Underbed light for safe entry and exit
- Massage mode for pain relief

SPECIFICATIONS

- Overall:
2050mm (L) x 1120mm (W) x 370mm (H)
- Platform Height: **370 - 690mm** - without mattress
- Product Weight: **114kg**
- Max. User Weight: **200kg**

WARRANTY

- Framework: **10 years**
- Electric Components: **3years**



ASK FOR THE BRAND

redgumbrand.com.au

Comfort Lux - RG66BS

Standard King Size Single Bed



Cordless Hand
Controller



Upholstery

FEATURES

- Multifunction - "Zero - G" and "Antisnore"
- Battery backup
- Underbed light for safe entry and exit
- Massage mode for pain relief
- 4 USB ports

SPECIFICATIONS

- Overall:
2020mm (L) x 1050mm (W) x 350mm (H)
- Height Adjustment: **350 - 450mm** - 5 settings
- Product Weight: **67kg**
- Max. User Weight: **200kg**

WARRANTY

- Framework: **10 years**
- Electric Components: **3years**



FUNCTIONS



Height Adjustable

Back Raise

Leg Raise

Zero Gravity

Anti Snore

ANTI-SNORE:

The adjustable bed head segment rises to an approximate 15 degree angle. This small movement can open up partially restricted airways.

ASK FOR THE BRAND

redgumbrand.com.au


Comfort Lux Mattress - RG66M200 & RG66M250

Medium and High Pressure Care Mattress



ICE COOL COVER

- Specifically designed for “hot” sleepers and those who prefer a cooler sleep experience.
- Manages wetness better through faster moisture wicking, airflow and ability to dry fast.
- The cooling UHMWPE fiber has an advanced moisture wicking structure that dries 2x faster than regular fabric.



Medium Pressure Care Mattress Specifications:

1050mm (W) x 2030mm (L) x 200mm (H) - Max. User Weight: 200kg

High Pressure Care Mattress Specifications:

1050mm (W) x 2030mm (L) x 250mm (H) - Max. User Weight: 200kg

Medium Pressure Care Mattress

- 20cm memory foam
- 3cm 30D gel infused memory foam
- 3cm 40D memory foam
- 14cm 25D support foam ventilated design
- Side wall: Three dimensional polypropylene with zipper in the middle
- Bottom: non-slip fabric

High Pressure Care Mattress

- 25cm memory foam
- 3cm 30D gel infused memory foam
- 5cm 40D memory foam
- 17cm 25D support foam ventilated design
- Side wall: Three dimensional polypropylene with zipper in the middle
- Bottom: non-slip fabric

Warranty: 5years (Against faulty workmanship and / or faulty materials)


ASK FOR THE BRAND

redgumbrand.com.au

Pressure Mapping Test on our Medium and High Pressure Care Mattress

It is an unfortunate fact that estimated pressure injury prevalence in Australian hospitals ranges from 9.5% to 17.6%, and in nursing home and long-term care setting is estimated to be around 8.9%. These wounds bring pain, associated risk for serious infection and increased healthcare utilization.

With our pressure mapping technology it is possible to objectively measure interface pressure and provide computer generated data that is **not based on opinion**. By means of our pressure mapping technology, it is possible to reduce interface pressure by selecting an appropriate pressure-relieving surface or use adjustable and active mattress to improve the quality of life.



Comfort Lux - Medium Pressure Care Mattress



Comfort Lux - High Pressure Care Mattress



Our pressure Care Mattresses can be used in circumstances where the user may be required to spend a large number of hours in bed per day (on top of when they are sleeping).

Research has shown, when comparing a regular mattress to a pressure care mattress:

- A significant decrease in pressure ulcer complications
- A decrease in time taken to fall asleep
- An increase in the amount of "Deep Sleep Time"

Pressure injuries can be caused by number of factors, including but not limited to:

- Impaired mobility
- Impaired activity
- Compromised skin integrity
- Compromised or reduced blood supply to pressure points

Surveys conducted throughout Hospital and Care Facilities in Australia have shown that 48.8% of pressure ulcers are found on the lower limbs (from the knee to toes) and 33.3% around the pelvis and lower back areas. A "Zero G Position" which is available on some High-Low beds is shown to significantly reduce the pressure ulcers caused in these areas of concern.

Our COMFORT LUX Pressure Care Mattress has been designed as a solution for users with Medium to High pressure care needs and to enable high levels of pressure relief, whilst still providing a high level of comfort in their living environment.

ASK FOR THE BRAND

redgumbrand.com.au

REDGUM BRAND
Daily Living and
Lifestyle Products



LUX COMFORT MAXIMUM FUNCTIONALITY

Please contact REDGUM Brand if you require further information on this product.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand: 1 Business way Malaga Western Australia
6090 Part of the For-de Group Pty Ltd

ASK FOR THE BRAND

redgumbrand.com.au