



INJECTION RED

CHARCOAL GREY

YELLOW



FEATURES

- Easy folding process for storage and transport
- Patented easy fold mid back seat
- Easy folding base
- Lifting handles makes picking up the scooter easy
- Multi position steering tiller
- Straight handle with finger operated controls
- Rear wheel drive
- Electromagnetic brake
- Flat free tyres
- Side reflectors

SPECIFICATIONS

Overall: 1010mm (L) x 445mm (W) x 890mm (H)

Folded: 470mm (L) x 445mm (W) x 775mm (H)

Turning Circle: 1040mm

Speed: 6km/h

Maximum range: Up to 18km

Safe climb angle: 6°

Maximum climb angle: 6°

Motor: 0.330Hp high torque motor (2 Pole)


Batteries: 2 x 10 Amp minimum (supplied)

Charger: 2 Amp

Tyres (Flat free): Ø155mm Front / Ø205mm Rear

Product Weight: 27kg - with batteries

Maximum User Weight: 90kg - on flat ground



0° - 90kg
6° - 75kg