

Jade Lightweight Steel Wheelchair

RG18SD / RG20SD



Features

- Lightweight steel construction
- Powder coated folding frame
- Padded arm rests
- Lift up armrests
- 24" removable solid rear tyres
- 8" Front rubber castors
- Hard wearing nylon upholstery
- Easy to operate swing away leg rests
- Manual brake on rear wheels
- Easy remove rear wheels
- Available in Metallic Green

Dimensions

- Seat Width 18" (450mm) RG18SD
- Seat Width 20" (510mm) RG20SD
- Seat Depth 430mm
- Seat Height 400mm
- Max User Capacity - 110kg
- Wheelchair Weight 18" - 13.0kg
- Wheelchair Weight 20" - 13.5kg

Options Available

- Anti Tip Wheels
- Elevating Leg rests (Chrome)

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@redgumbrand.com.au