

COUGAR 10 LIFT

BURGUNDY RED

G424ML

FEATURES

- High back contoured seat with headrest
- Electric power lift height / angle adjustable seat
- Lift up & adjustable width / height arm rests
- Swing away legrest
- Mid wheel drive
- Finger control operation
- High / Low speed settings
- Battery power indicator
- Horn
- Ideal for indoor and outdoor use on even ground
- Tie down points for transport

- Attachment options available:
 Elevating swing away legrests
 Stump support

SPECIFICATIONS

Overall: 920mm (L) x 620mm (W) x 1240mm (H)

Seat height range: 580 - 910mm

Turning circle: 950mm

Speed: 7 - 8km/h

Maximum range: Up to 40km

Safe climb angle: 9°

Maximum climb angle: 9°

Motor: 2 x 0.8Hp high torque motors (2 Pole)


Batteries: 2 x 50 Amp minimum (supplied)

Charger: 5 Amp

Tyres: Front: Ø150mm / Middle: Ø260mm / Rear: Ø125mm

Product Weight: 114kg - with batteries

Maximum User Weight: 136kg - on flat ground



0° - 136kg
 6° - 113kg
 9° - 91kg