

REDGUM BRAND
Daily Living and
Lifestyle Products

Deluxe Mobile Tilting Comfort Chair

RGCCAD



Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@redgumbrand.com.au

REDGUM BRAND
Daily Living and
Lifestyle Products

REDGUM Brand is a 100% Australian Owned & Operated

Deluxe Mobile Tilting Comfort Chair

RGCCAD



Features

- Memory foam seat base
- Height adjustable extra long armrests
- Armrests adjust for easy transfer
- Seat base is adjustable to 12 degrees
- Seat backrest is adjustable to 30 degrees
- Gas assisted seat adjustment operation
- Dual controls for patient & attendant
- Attendant operated tilt & space function
- Removable adjustable headrest
- Side supports that swing away
- Centralised foot operated brake system
- Retractable foot plate with damper
- Removable meal table as standard

Dimensions

- Seat width 500mm
- Seat depth 490mm
- Seat height 550mm
- Seat base angle -8 to 12°
- Back rest height 780mm
- Back rest width 470mm
- Back rest angle to 30°
- Armrest height 545 to 800mm
- Overall width 660mm
- Overall height 1320mm
- Overall length 900mm
- Max user weight 140kg
- Weight of chair 35.0kg

Optional Extras

- I V Pole
- Document holder
- Urine bag holder
- Oxygen bottle holder



Removable meal table as standard