

**REDGUM** BRAND  
Daily Living and  
Lifestyle Products

---

## Deluxe Comfort Chair

RGCCN



---

Phone: 08 9248 4180  
Fax: 08 9248 9452  
E-mail: [orders@redgumbrand.com.au](mailto:orders@redgumbrand.com.au)

**REDGUM** BRAND  
Daily Living and  
Lifestyle Products

*REDGUM Brand is a 100% Australian Owned & Operated*

## Deluxe Comfort Chair

RGCCN



### Features

- Memory foam seat base
- Contoured seat design
- Integrated leg rest when reclined
- Gas assisted seat adjustment operation
- Rear attendant handle for easy movement
- Height adjustable extra long armrests
- Armrests adjust for easy transfer
- Seat backrest adjustable 10 to 45 degrees
- Removable headrest bolster
- Locking castor on each leg
- Durable vinyl seat covering
- Meal table comes as standard

### Dimensions

- Seat width 535mm
- Seat depth 430mm
- Seat height 540mm
- Seat base angle 10°
- Back rest height 780mm
- Back rest width 470mm
- Back rest angle to 45
- Armrest height 535 to 760mm
- Overall width 670mm
- Overall height 1225mm
- Overall length 840mm
- Max user weight 130kg
- Weight of chair 28.1kg