

# Shoprider

## 778ER

ASK FOR THE BRAND SHOPRIDER



**SHOPRIDER**<sup>®</sup>  
*mobility*

[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

## 778ER

**SHOPRIDER**<sup>®</sup>  
www.shoprider.com.au *mobility*

The SHOPRIDER™ 778ER is at home both indoors and outside. It is a highly manoeuvrable scooter suited to tight areas within the house or out shopping. The 778ER has great features as standard including angled foot platform for extra leg room, fully rotating lockable high back captains seat, adjustable steering tiller and fold up arm rests. The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

### SPECIFICATIONS

### 778ER

Length	1120mm
Width	600mm
Height	1160mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	0.80Hp
Maximum Speed	8Kph
Battery Capacity (Minimum Amp)	2 x 35Amp
Maximum Range #	up to 26km
Wheel Size	260mm
Turning Radius	1000mm
Weight (w/o batteries)	79kg

### FEATURES

High back Captain's Seat  
 Loop Steering Tiller  
 Adjustable Steering Tiller  
 Adjustable Width Arm Rests  
 Sliding Rotating Seat  
 Variable Seat Height  
 Battery Power indicator  
 Head Lights  
 Rear Lights & indicators  
 Hazard Lights  
 Shopping Basket  
**Comprehensive Warranty**


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

**WARNING:** Do not exceed 10kph on public roads or paths

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



**SHOPRIDER**<sup>®</sup>  
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
 Or visit our Web Site at www.shoprider.com.au