

Shoprider

778HD

ASK FOR THE BRAND SHOPRIDER



SHOPRIDER[®]
mobility
www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

778HD

SHOPRIDER[®]
www.shoprider.com.au *mobility*

The SHOPRIDER™ 778HD is another fantastic scooter from the extensive range of scooters and powerchairs by SHOPRIDER. The 778HD is the workhorse of the three wheel range of scooters offered it has all the features that expect from Shoprider, fully rotating Airflow comfort seat,, suspension, adjustable steering tiller and the added bonus of flat free tyres as standard. This scooter also a powerful 2hp motor. The Reliability of this impressive scooter by SHOPRIDER is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

Length	1230mm
Width	640mm
Height	1250mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	182Kgs
Charger	4Amp
Motor Size	2.0Hp
Maximum Speed	12Kph
Battery Capacity (Minimum Amp)	2 x 50Amp
Maximum Range #	up to 30km
Wheel Size	260/330mm
Turning Radius	1290mm
Weight (w/o batteries)	76kg

FEATURES

High Back Seat
Adjustable Steering Tiller
Adjustable Width Arm Rests
Fold Arm Rests
Sliding Rotating Seats
Front Shopping Basket
Variable Seat height
Horn
Battery Level Warning Light
Flat Free Tyres
Front & Rear Suspension
Front Lights
Anti Tips Wheels
Comprehensive Warranty


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING: Do not exceed 10kph on public roads or paths

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



SHOPRIDER[®]
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au