

# Shoprider

## 888SE

ASK FOR THE BRAND SHOPRIDER



Available Colours



**SHOPRIDER**<sup>®</sup>  
*mobility*

[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

## 888SE

www.shoprider.com.au *mobility*

The SHOPRIDER™ 888SE is packed with all the features of the larger models with the advantage of its compact size. The 888SE equipped with large batteries and has the following features, 260mm tyres for safety & comfort, fully rotating lockable seat, adjustable steering tiller and fold up arm rests. The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

### SPECIFICATIONS

### 888SE

Length	1280mm
Width	600mm
Height	1070mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	8Amp
Motor Size	0.75Hp
Maximum Speed	8Kph
Battery Capacity (Minimum Amp)	2 x 50Amp
Maximum Range #	up to 40km
Wheel Size	260mm
Turning Radius	1500mm
Weight (w/o batteries)	76kg

### FEATURES

High Back Seat  
 Adjustable Steering Tiller  
 Adjustable Width Arm Rests  
 Fold Arm Rests  
 Sliding Rotating Seats  
 Front Shopping Basket  
 Variable Seat height  
 Horn  
 Battery Level Warning Light  
 Front & Rear Bumper Pads  
 Front & Rear Suspension  
 Front Lights  
 Anti Tips Wheels  
**Comprehensive Warranty**

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia and is manufactured to comply with all relevant international ISO standards.


**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

**WARNING:** Do not exceed 10kph on public roads or paths

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

### 888SE



Your Authorised dealer

For your nearest Distributor call 1800 672 633  
 Or visit our Web Site at www.shoprider.com.au