

Shoprider COUGAR

ASK FOR THE BRAND SHOPRIDER



COUGAR 14



COUGAR 10

SHOPRIDER[®]
mobility
www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER®

www.shoprider.com.au

mobility

COUGAR

The SHOPRIDER™ COUGAR is ergonomically engineered to provide the ultimate combination of style and comfort. **Cougar** with its powerful smooth motors, mid wheel maneuverability, and comfortable seating provides a safe and comfortable ride. The smooth operation combined with mid wheel drive means easy operation under various conditions. It's ideal for both inside and outside use on even ground and can travel long distances between charges.

SPECIFICATIONS	COUGAR 10	COUGAR 14
Length	920mm	920mm
Width	620mm	690mm
Height	1120mm	1180mm
Battery Capacity (Minimum Amp)	2 x 40Amp	2 x 40Amp
#Maximum Range	up to 30km	up to 40km
Maximum Load Single User*	136Kgs	136Kgs
Charger	4Amp	8Amp
Motor Size	0.8 Hp	1.34 Hp
Maximum Speed	8Kph	9.6Kph
*Safe Climbing Angle	8 Degrees	8 Degrees
*Max Climbing Angle	9 Degrees	12 Degrees
Turning Radius	570mm	600mm
Wheel Size:		
- Front Caster	2 x 125mm	2 x 125mm
- Middle (drive)	2 x 260mm	2 x 360mm
- Rear caster	2 x 100mm	2 x 100mm

FEATURES:

High Back Captain's Seat
Adjustable Angle and Height Seat
Adjustable Armrests
Finger Control Operation
Tie Down Points for Transport
Swing Away Elevating Available
Legs Rest Available
Stump Supports Available


Driver weight may exceed weight of unit. Reduce speed when turning

IMPORTANT

Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



SHOPRIDER®
www.shoprider.com.au

mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au