



www.shoprider.com.au

mobility

THE GECKO RANGE OF PORTABLE MOBILITY SCOOTERS

The Gecko SL-3 & SL-4 are cost effective mobility scooters with a modern sleek design, comfortable seat and the portability to transport easily. The Gecko has a comfortable seat with a built in storage pocket. The reliability of this impressive little scooter is further enhanced by the rear wheel drive and electromagnetic brake.



SL-3 Specifications

Overall Dimensions:	L x 950 W x 510 H x 920
Total Weight:	28 kgs
Weight of heaviest part:	11 kgs
Battery pack weight (10amp)	7 kgs
Charger (off board)	2Amps
Max single user weight (on flat road)	105 kgs
Maximum speed	6 Km/Hr
Range on flat ground	12 Km
Turning Radius	835mm

Specifications subject to change without notice



SL-4 Specifications

Overall Dimensions:	L x 970 W x 510 H x 920
Total Weight:	31 kgs
Weight of heaviest part:	11 kgs
Battery pack weight (10amp)	7 kgs
Charger (off board)	2Amps
Max single user weight (on flat road)	110 kgs
Maximum speed	6 Km/Hr
Range on flat ground	12 Km
Turning Radius	1080mm

Specifications subject to change without notice

FEATURES: Light weight design, Comprehensive warranty, Solid tyres, Anti tip wheels, Battery meter, Horn, Easy to use controls, Drop down armrests, Available in Red or Blue

For your nearest distributor call
1800 672 633 or visit our website at
www.shoprider.com.au

Your Authorised Dealer: