

# Shoprider

## LITTLE BEAUTY

ASK FOR THE BRAND SHOPRIDER



Available Colours



**SHOPRIDER**<sup>®</sup>  
*mobility*  
[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

**SHOPRIDER**<sup>®</sup>  
www.shoprider.com.au *mobility*

## LITTLE BEAUTY

The SHOPRIDER™ **LITTLE BEAUTY** is a versatile portable scooter that offers the performance and features that you expect from a SHOPRIDER™ scooter. This portable scooter is easily disassembled for transport with NO requirement for tools. The **LITTLE BEAUTY** has been designed to offer a number of improved features like LED battery power display, larger leg room, upgradeable battery box for larger batteries and comfortable swivel seat. The reliability of this impressive portable scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

### SPECIFICATIONS

### LITTLE BEAUTY

### FEATURES

Length	1110mm
Width	540mm
Height	925mm
*Safe Climbing Angle	6 Degrees
*Max Climbing Angle	6 Degrees
Maximum Load Single User*	113Kgs
Charger	2 Amp
Motor Size	0.3Hp
Maximum Speed	6Kph
Battery Capacity (Minimum Amp)	12-21Amp
Maximum Range #	up to 23km
Wheel Size	200mm
Turning Radius	1170mm
Weight (with battery)	43kg (12 amp Batts)

Mid Back Seat  
Adjustable Steering Tiller  
Adjustable Width Arm Rests  
Rotating Seat  
Folding Arm Rests  
Variable Seat height  
Front Shopping Basket  
Horn  
LED Battery Level indicator  
Upgradeable Battery Size  
Anti Tips Wheels  
**Comprehensive Warranty**

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia and is manufactured to comply with all relevant international ISO standards.


**WARNING:** Do not exceed 10kph on public roads or paths.

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

### LITTLE BEAUTY



**SHOPRIDER**<sup>®</sup>  
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
Or visit our Web Site at www.shoprider.com.au