

Shoprider

Rocky 4

ASK FOR THE BRAND SHOPRIDER



ROCKY 4

SHOPRIDER[®]

www.shoprider.com.au

mobility

Available Colours



RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER
www.shoprider.com.au *mobility*

Rocky 4

The SHOPRIDER™ **Rocky 4** is an outstanding model from the extensive range of mobility scooters available from SHOPRIDER. The **Rocky 4** must be seen to be believed, it redefines the concept of reliable mobility. The **Rocky 4** comes with 330mm wheels, Airflow Comfort Seat and all the features you expect from the SHOPRIDER range, fully rotating lockable seat, suspension, adjustable steering tiller and multi position lift up armrests.

The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS	ROCKY 4
Length	1430mm
Width	680mm
Height	1300mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
*Maximum Load <u>Single user</u>	225Kgs
Charger	8Amp
Motor Size	2.0Hp
Maximum Speed	12Kph
Battery capacity min	2 x 60 Amp
Maximum range#	up to 50km
Wheel Size	330mm
Turning Radius	1800 mm
Weight (with battery)	146kg

FEATURES:


- Airflow Comfort Seat
- Loop Steering Tiller
- Gas Adjustable Steering Tiller
- Adjustable width arm rests
- Sliding Rotating Seat
- Variable Seat Height
- Battery Power indicator
- Double Head Lights
- Rear Lights & Indicators
- Hazard Lights
- Front & Rear Bump Pads
- High & Low Speed
- Shopping Basket
- Comprehensive Warranty**

WARNING: Do not exceed 10kph on public roads or paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)



Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au