

# SHOPRIDER®

Mobility

PRODUCT UPDATE !



## ROCKY 6 FAT TYRE ELECTRIC SCOOTER

889XLSBNB - RED

889XLSBNB - BLUE

### FEATURES

#### - **NEW LED HEADLIGHTS!**

- Double pillow comfort high back seat with headrest
- Gas adjustable delta steering tiller
- Loop handle with finger operated controls
- High / low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + rear (adjustable) suspension
- Twin Headlights (LED)
- Rear lights
- Indicators
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up adjustable width arm rests
- Removable shopping basket
- Front + rear bump pads

### SPECIFICATIONS (All dimensions shown are nominal)

Overall: **1430mm (L) x 680mm (W) x 1300mm (H)**  
Turning Circle: **1800mm**  
Speed: **10kmph**  
Maximum range: **Up to 50km**  
Safe climb angle: **9°**  
Maximum climb angle: **12°**  
Motor: **2.0Hp high torque motor**  
Batteries: **2 x 60 amp minimum (2 x 75 amp supplied)**  
Charger: **8 amp**  
Tyres: **Ø330mm pneumatic - Front + Rear**  
Product Weight: **147kg - with 60 amp batteries**  
Maximum User Weight: **225kg - on flat ground**



LED Headlights

#### IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

Please contact **SHOPRIDER** if you require further information on this product.

**P: 08 9248 4180**

**F: 08 9248 9458**

**E: [orders@shoprider.com.au](mailto:orders@shoprider.com.au)**

SHOPRIDER - 1 Business way, Malaga, Western Australia 6090  
Part of the **FOR-DE Group** Pty Ltd.

**[shoprider.com.au](http://shoprider.com.au)**