

SHOPRIDER®

Mobility

PRODUCT UPDATE !



LED Headlights



SHOPRIDER strives to provide the highest quality products possible that are subject to change without notice.

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

ALLROUNDER ELECTRIC SCOOTER

9LSP - RED

9LSP - BLUE

9LSP - CHAMPAGNE

FEATURES

- **NEW LED HEADLIGHTS!**

- Comfort high back seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + rear (adjustable) suspension
- Twin Headlights (LED)
- Rear lights
- Indicators
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up adjustable width arm rests
- Removable shopping basket
- Front chrome bumper

SPECIFICATIONS (All dimensions shown are nominal)

Overall: **1300mm (L) x 640mm (W) x 1200mm (H)**

Turning Circle: **1440mm**

Speed: **10kmph**

Maximum range: **Up to 40km**

Safe climb angle: **9°**

Maximum climb angle: **12°**

Motor: **1.5Hp high torque motor**

Batteries: **2 x 40 amp minimum**

Charger: **8 amp**

Tyres: **Ø330mm pneumatic - Front + Rear**

Product Weight: **99kg - with 40 amp batteries**

Maximum User Weight: **159kg - on flat ground**

Please contact **SHOPRIDER** if you require further information on this product.

P: 08 9248 4180

F: 08 9248 9458

E: orders@shoprider.com.au

SHOPRIDER - 1 Business way, Malaga, Western Australia 6090
Part of the **FOR-DE Group** Pty Ltd.

shoprider.com.au