

# Shoprider

## Little Ripper

ASK FOR THE BRAND SHOPRIDER



**GK9**

**SHOPRIDER**<sup>®</sup>  
*mobility*  
[www.shoprider.com.au](http://www.shoprider.com.au)

Available Colours



RELIABLE MOBILITY GUARANTEED

# Shoprider

## Little Ripper



The SHOPRIDER™ **Little Ripper** is a fully portable model that offers exceptional functionality for each user. The Little Ripper has many features found on the larger scooters with the added benefit of its easy pull apart design. Features include puncture proof tyres, mid back fully rotating seat, adjustable steering tiller, multi position lift up armrests & 21amp batteries. The little Ripper is suitable for indoor and outdoor usage and is ideal for easy transport when pulled apart. The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

### SPECIFICATIONS

SPECIFICATIONS	GK9
Length	1050mm
Width	560mm
Height	880mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	0.8Hp
Maximum Speed	6kmph
Battery Capacity (Minimum Amp)	21 x 21 Amp
Maximum Range #	up to 20km
Wheel Size	225mm
Turning Radius	1150mm
Weight (with battery)	56kg
Weight (Battery Each)	7.5kg (15kg Total)
Weight (Chair)	10.3kg
Weight (Rear Transaxle)	14.3kg
Weight (Main Chassis)	16.4kg

### FEATURES

- Mid Back Seat
- T Bar Steering Tiller
- Adjustable Steering Tiller
- Adjustable width arm rests
- Sliding Seat
- Variable Seat Height
- Battery Power indicator
- Front Bump Bar
- Pull Apart Design for Transport
- Shopping Basket
- Comprehensive Warranty**


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

**WARNING:** Do not exceed 10kph on public roads or paths.

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



Your Authorised dealer

For your nearest Distributor call 1800 672 633  
Or visit our Web Site at [www.shoprider.com.au](http://www.shoprider.com.au)