

Shoprider Getabout

ASK FOR THE BRAND SHOPRIDER



GETABOUT

SHOPRIDER[®]
mobility
www.shoprider.com.au

Available Colours



RELIABLE MOBILITY GUARANTEED

Shoprider

Getabout

SHOPRIDER[®]
www.shoprider.com.au *mobility*

The SHOPRIDER™ **Getabout** has been designed to be easily folded and transported, ideal for people wanting to travel via vehicle, aeroplane or ship. The **Getabout** features a patented easy fold seat mechanism to make the folding process very quick and when you combine this with the travel handle and unidirectional castors you have a great scooter for your portable needs. The **Getabout** is another innovative scooter from Shoprider.

The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS	GETABOUT
Length	1010mm (470mm)
Width	445mm
Height	890mm
*Safe Climbing Angle	6 Degrees
*Max Climbing Angle	6 Degrees
*Maximum Load <u>Single user</u>	90Kgs
Charger	2Amp
Motor Size	0.330Hp
Maximum Speed	6Kph
Battery capacity min	2 x 10 Amp
Maximum range#	up to 18km
Wheel Size	155/205mm
Turning Radius	915 mm
Weight (with battery)	27kg

FEATURES:


- Patented Easy Fold Seat Mechanism
- Easy Folding Base
- Multi Position Tiller
- Flat Free Tyres
- Travel Handle
- Unidirectional Castors
- 10amp Batteries
- Comprehensive Warranty**

WARNING: Do not exceed 10kph on public roads or paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



*(Single user)

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

SHOPRIDER[®]
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au