

COUGAR 14

BURGUNDY RED


G424L

FEATURES

- High back contoured seat with headrest
 - Height / angle adjustable seat
 - Lift up & adjustable width / height arm rests
 - Fold up adjustable height / depth footrest
 - Mid wheel drive
 - Finger control operation
 - High / Low speed settings
 - Battery power indicator
 - Horn
 - Ideal for indoor and outdoor use on even ground
 - Tie down points for transport
- Attachment options available:
 Swing away legrests
 Elevating swing away legrests
 Stump support

SPECIFICATIONS

- Overall: 990mm (L) x 690mm (W) x 1180mm (H)
 Turning Circle: 1020mm
 Speed: 9.6km/h
 Maximum range: Up to 40km
 Safe climb angle: 8°
 Maximum climb angle: 12°
 Motor: 2 x 1.34Hp high torque motors (2 Pole)
 Batteries: 2 x 50 Amp minimum (supplied)
 Charger: 8 Amp
 Tyres: Front: Ø150mm / Middle: Ø360mm / Rear: Ø125mm
 Product Weight: 118kg - with batteries
 Maximum User Weight: 136kg - on flat ground



0° - 136kg
 6° - 113kg
 9° - 91kg
 12° - 75kg