

BURGUNDY RED



UL8-W

FEATURES

- Powerchair pulls apart for storage and transport
- Comfortable mid back contoured seat
- Height / angle adjustable seat
- Lift up & adjustable width / height arm rests
- Fold up adjustable height / depth footrest
- Mid wheel drive
- Finger control operation
- High / Low speed settings
- Battery power indicator
- Horn
- Ideal for indoor and outdoor use on even ground
- Safe to transport on planes

SPECIFICATIONS

Overall: 725mm (L) x 545mm (W) x 1180mm (H)
Turning Circle: 755mm
Speed: 6km/h
Maximum range: Up to 12km
Safe climb angle: 10°
Maximum climb angle: 10°
Motor: 2 x 0.5Hp high torque motors (2 Pole)
Batteries: 10 Amp minimum (Lithium Iron Phosphate)
Charger: 2 Amp
Tyres: Front: Ø155mm / Rear: Ø200mm
Maximum User Weight: 113kg - on flat ground

Product Weight: 41kg - with batteries
Seat with Controller Weight: 15.7kg
Front Chassis with Battery Weight: 12.7kg
Rear Chassis with Drive Train Weight: 12.6kg



0° - 113kg
6° - 100kg
10° - 90kg