

BURGUNDY RED



FEATURES

778HD

- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear suspension
- Headlights
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

SPECIFICATIONS

Overall: 1230mm (L) x 640mm (W) x 1250mm (H)

Turning Circle: 1250mm

Speed: 10km/h

Maximum range: Up to 40km

Safe climb angle: 9°

Maximum climb angle: 9°

Motor: 2.0Hp high torque motor (2 Pole)

Batteries: 2 x 50 Amp minimum (supplied)

Charger: 5 Amp

Tyres (Pneumatic): Ø260mm Front / Ø330mm Rear

Product Weight: 102kg - with batteries

Maximum User Weight: 182kg - on flat ground

