

BURGUNDY RED



778SL

FEATURES

- High back contoured seat with headrest
- Gas adjustable steering tiller with modesty panel
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear suspension
- Twin Headlights
- Rear lights
- Indicators
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

SPECIFICATIONS

Overall: 1230mm (L) x 590mm (W) x 1250mm (H)

Turning Circle: 1250mm

Speed: 10km/h

Maximum range: Up to 35km

Safe climb angle: 9°

Maximum climb angle: 9°

Motor: 1.2Hp high torque motor (2 Pole)


Batteries: 2 x 50 Amp minimum (supplied)

Charger: 5 Amp

Tyres (Pneumatic): Ø260mm Front / Ø330mm Rear

Product Weight: 88kg - with batteries

Maximum User Weight: 136kg - on flat ground



0° - 136kg
6° - 113kg
9° - 91kg