


BURGUNDY RED



0° - 225kg
6° - 180kg
9° - 159kg
12° - 136kg

889D

FEATURES

- Comfortable large wide bench seat - 800mm
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Twin Headlights (LED)
- Rear lights
- Indicators
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding seat
- Lift up & adjustable width arm rests
- Removable shopping basket

SPECIFICATIONS

Overall: 1480mm (L) x 910mm (W) x 980mm (H)
 Turning Circle: 2550mm
 Speed: 10km/h
 Maximum range: Up to 40km
 Safe climb angle: 9°
 Maximum climb angle: 12°
 Motor: 2.0Hp high torque motor (2 Pole)
 Batteries: 2 x 75 Amp minimum (supplied)
 Charger: 8 Amp
 Tyres (Pneumatic): Ø330mm - Front + Rear
 Product Weight: 159kg - with batteries
 Maximum User Weight: 225kg - on flat ground
SINGLE PERSON USE ONLY