



9LSP

FEATURES

- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- LED Headlight
- Rear lights
- Indicators
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

SPECIFICATIONS

Overall: 1300mm (L) x 640mm (W) x 1200mm (H)
 Turning Circle: 1440mm
 Speed: 10km/h
 Maximum range: Up to 40km
 Safe climb angle: 9°
 Maximum climb angle: 12°
 Motor: 1.5Hp high torque motor (4 pole)
 Batteries: 2 x 50 Amp minimum (supplied)
 Charger: 8 Amp
 Tyres (Pneumatic): Ø330m - Front + Rear
 Product Weight: 99kg - with batteries
 Maximum User Weight: 159kg - on flat ground

