

LITTLE BEAUTY

BURGUNDY RED

BLUE



QT8

FEATURES

- Scooter pulls apart for storage and transport
- Comfortable mid back contoured seat
- Fold down adjustable steering tiller
- Straight handle with finger operated controls
- Rear wheel drive
- Electromagnetic brake
- Battery power indicator
- Upgradable battery size
- Horn
- Removable rotating / height adjustable seat
- Lift up & adjustable width arm rests
- Lifting handles makes picking up the scooter easy
- Removable shopping basket

SPECIFICATIONS

Overall: 1110mm (L) x 540mm (W) x 925mm (H)
 Turning Circle: 1130mm
 Speed: 6km/h
 Maximum range: Up to 18km
 Safe climb angle: 6°
 Maximum climb angle: 6°
 Motor: 0.3Hp high torque motor (2 Pole)
 Batteries: 2 x 12-21 Amp minimum (12 Amp supplied)
 Charger: 2 Amp
 Tyres (Flat free): Ø200mm - Front + Rear
 Maximum User Weight: 113kg - on flat ground

Product Weight: 45kg - with 12 amp batteries
 Battery Weight: 5.6kg each
 Seat Weight: 8.3kg
 Rear Transaxle Weight: 9.9kg
 Main Chassis Weight: 15.6kg

