

Shoprider

— ENDEAVOUR

ASK FOR THE BRAND SHOPRIDER



Available Colours



SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER[®]
www.shoprider.com.au *mobility*

ENDEAVOUR

The SHOPRIDER™ ENDEAVOUR is the latest scooter from SHOPRIDER™ and it is loaded with great features. The ENDEAVOUR has the all the standard SHOPRIDER™ features 330mm tyres for safety & comfort, fully rotating lockable seat, adjustable steering tiller, fold down arm rests plus powerful LED lights, usb outlet and modern new design body covers. The reliability of these impressive scooters by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

ENDEAVOUR

FEATURES

Length	1370mm
Width	630mm
Height	1360mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	159Kgs
Charger	8Amp
Motor Size - 4Pole	0.9Hp
Maximum Speed	10Kph
Battery Capacity (Minimum Amp)	40Amp
Maximum Range #	up to 32km
Wheel Size	330mm
Turning Radius	1750mm
Weight (with battery)	101kg

Large Captain's Seat
Adjustable Steering Tiller
Adjustable Width Arm Rests
Sliding Rotating Seat
Front Shopping Basket
Horn
USB power outlet
LED Battery Level indicator
LED Front Headlights
Front & Rear Suspension
Front & Rear Lights
Hazard Lights /indicators
High / Low Speed switch
Anti Tips Wheels
Comprehensive Warranty

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia and is manufactured to comply with all relevant international ISO standards.


WARNING: Do not exceed 10kph on public roads or paths.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

ENDEAVOUR



SHOPRIDER[®]
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au