

FPC



Presenting **Shoprider's** newest folding power chair the **FPC**. The **FPC** is packed with pioneering technology and style. The steel red frame, puncture proof tires and height adjustable arms are all standard with the **FPC**. The standard seat cushion on the **FPC** makes it a pleasure to drive.

Available Colors




Red

VSI Controller

Adjustable
Height Arm Rest

Footrest with
Length Adjustment

13" Wheel



Flip-Up Arm Rest




Adjustable Height Arm Rest



Adjustable Length
VSI Controller

**The FPC easily
folds to keep
you on the go!**



FPC

- Seat cushion for supreme comfort
- Available in 18" and 20" seat base
- Easy access on and off chair
- Adjustable length VSI controller
- Easy to disassemble for transportation

Accessories



Oxygen Bottle Holder



Elevating Leg Rest

FPC Specifications

Overall Dimensions (L' x W x H)	mm	900 x 620 x 1250
Number, Size of Tires (front; rear)	ea, in	2, 8; 2, 13
Weight Capacity	kg	Max. 136kg
Weight of Heaviest Part	kg	16.6kg
Total Unit Weight	kg	43kg
Battery Module Weight	kg	11.3kg each x 2
Battery Module Capacity	_V_Ah x Pcs	12V36Ah x 2
Charger		Off-board (4A)
Output Power of Motor	hp	0.6
Turning Radius	mm	650mm
Ground Clearance	mm	720
Maximum Speed	kph	6

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below.

a) The weight of occupant
b) Ground surface
c) Battery capacity and conditions
d) Type of charger
e) Ambient temperature
f) The way of driving
g) If the battery and mechanical moving parts fully break in
h) Etc.

SHOPRIDER

www.shoprider.com.au