

Shoprider COUGAR

ASK FOR THE BRAND SHOPRIDER



COUGAR POWER TILT

SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER[®]
www.shoprider.com.au *mobility*

COUGAR

The SHOPRIDER™ **COUGAR** power tilt has all the features that come standard on the **Cougar** range of powerchair. The power tilt model features a tilting seat function which is easily operated via the joy stick controller. giving the operator a maximum recline of 45 degrees for greater user comfort. This impressive powerchair will provide smooth operation with its mid wheel drive and powerful motors. It's ideal for both inside and outside use on even ground and can travel long distances between charges.

SPECIFICATIONS

G424MT

FEATURES:

Length	1090mm
Width	620mm
Height	1240mm
Battery Capacity (Minimum Amp)	2 x 50Amp
#Maximum Range	up to 40km
Maximum Load Single User*	136Kgs
Charger on board	5Amp
Motor Size	0.80 Hp
Maximum Speed	7-8Kph
*Max Climbing Angle	9 Degrees
Turning Radius	600mm
Wheel Size:	
- Front Caster	2 x 125mm
- Middle (drive)	2 x 250mm
- Rear caster	2 x 100mm
Maximum Seat Recline	45 degree

Large Captain's Seat
Electric power tilt seat
Adjustable arm rests
Swing away leg rest
Sleek body design
Color change panels
Smooth operation
Free wheel device
Variable seat height
Controller interchangeable both arms
Off Board charger


Driver weight may exceed weight of unit. Reduce speed when turning

IMPORTANT

Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



SHOPRIDER[®]
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au