

# Shoprider VENICE

ASK FOR THE BRAND **SHOPRIDER**



**SHOPRIDER**<sup>®</sup>  
*mobility*

[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

## VENICE

www.shoprider.com.au

mobility

The SHOPRIDER™ VENICE powerchair features World Leading Battery Technology to reduce the overall weight of the unit. The Venice is a versatile powerchair ideally suited to the needs of the portable user. The Venice can be dismantled for transport by using the quick disconnect lever in seconds. The Venice features a comfortable mid back seat, adjustable height armrests and when combined with the powerful motors it will deliver outstanding safety & comfort.

### SPECIFICATIONS

### VENICE

### FEATURES:

Length	725mm	Lithium Iron Phosphate Battery
Width	545mm	Mid Back Captain Seat
Height	1180mm	Colour Change Body Covers
Battery Capacity -Lithium Iron Phosphate	10 Amp	Variable seat height
#Maximum Range	up to 16km	Adjustable arm rests
Maximum Load Single User*	113Kgs	Controller interchangeable
Charger on board	2Amp	Free wheel device
Motor Size	500 watt	Easy Dismantle for Transport
Maximum Speed	6Kph	
*Safe Climbing Angle	10 Degrees	
*Max Climbing Angle	10 Degrees	
Turning Radius	400mm	
Wheel Size:		
- Front Caster	2 x 155mm	
- Rear Caster	2 x 200mm	

Driver weight may exceed weight of unit. Reduce speed when turning

### IMPORTANT

Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



10'(17)-90kg  
6'(10)-100kg  
0'(0)- 113kg

\*(Single user)

**SHOPRIDER®**  
www.shoprider.com.au

mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
Or visit our Web Site at www.shoprider.com.au