

**SHOPRIDER®**  
*mobility*


**VIKING**  
779XLS



INJECTION RED

Ask for the brand **Shoprider!**

INJECTION RED



0° - 185kg  
6° - 136kg  
9° - 91kg



## FEATURES

779XLS

- High back contoured seat with headrest
- Sliding / Rotating / Height adjustable seat
- Lift up + tilt adjustable arm rests
- "Slow when turning" safety feature
- Elevated handlebar with R/H paddle control
- Full function LCD display
- High / Medium / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front disc handbrake
- Front + Rear suspension
- LED Headlight - low / high beam
- Rear activated brake lights
- Indicators
- Hazard lights
- Large 75 amp battery capacity
- Battery power indicator

## SPECIFICATIONS

Overall: 1550mm (L) x 700mm (W) x 1305mm (H)

Turning Circle: 1650mm

Speed: 10km/h

Maximum range: Up to 54km

Safe climb angle: 9°

Maximum climb angle: 9°

Motor: 2.0Hp high torque motor (4 Pole)

Batteries: 2 x 75 Amp minimum (supplied)

Charger: 8 Amp

Tyres (Pneumatic): Ø406mm Front / Rear (3.00-10)

Product Weight: 140kg - with batteries

Maximum User Weight: 185kg - on flat ground