

**SHOPRIDER**<sup>®</sup>  
mobility

**VIKING**



0° - 185kg  
6° - 136kg  
9° - 91kg

## FEATURES

**779XLS**

- High back contoured seat with headrest
- Sliding / rotating / height adjustable seat
- Lift up + tilt adjustable arm rests
- "Slow when turning" safety feature
- Elevated handlebar with R/H paddle control
- NEW Full function, colour LCD HD display
- High / medium / low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front disc handbrake
- Front + rear suspension
- LED Headlight - low / high beam
- Daytime running light
- Rear LED brake lights
- Front and rear indicators
- Hazard lights
- Large 75 amp battery capacity
- Battery power indicator

Rear fenders available in  
*Storm Grey* or *Pearl White*



**NEW DESIGN**


P: **08 9248 4180**

E: [orders@for-de.com.au](mailto:orders@for-de.com.au)


**ASK FOR THE BRAND**

[shoprider.com.au](http://shoprider.com.au)


IMPORTANT NOTE: It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.




Overall Length:  
**1550mm (L)**




Seat Width:  
**515mm**



Backrest Height:  
**465mm**




Armrest Length:  
**335mm**




Headrest Height  
(from seat): **4 settings**


**710mm**  
**690mm**  
**660mm**  
**640mm**




Overall Width:  
**700mm (W)**




Seat Depth:  
**510mm**




Backrest Width:  
**480mm**




Armrest Width:  
**60mm**




Headrest Width:  
**250mm**




Overall Height  
(Lowest seat setting):  
**1305mm (H)**




Seat Rotation:  
**90° Left & Right**



Backrest Angle:  
**Folded down to 109°**




Distance between  
armrests: **485mm**




Seat edge to centre of tiller:  
(adjustable) **425 - 290mm**

**ASK FOR THE BRAND**


**shoprider.com.au**




Seat Height **720mm**  
(from floor): **695mm**  
**3 settings** **670mm**




Ground Clearance: **150mm**




180° Turning Circle: **3120mm**




Anti-tip Wheel  
Ground Clearance: **30mm**




Safe & Maximum  
climb angle: **9°**



Seat Height **530mm**  
(from footwell): **505mm**  
**3 settings** **480mm**




Speed: **10km/h**




Tyres (Pneumatic):  
**Front Ø406mm**  
**Rear (3.00-16)**




Max. User Weight:  
**185kg - on flat ground**



Space Between Handlebars:  
**550mm**




Charger:  
**10 Amp**



Maximum range:  
**Up to 50km**




Batteries:  
**2 x 75 Amp**  
**min. (supplied)**



Motor:  
**2.0Hp high torque**  
**motor (4 Pole)**



Product Weight:  
**140kg - with batteries**



Screen Size:  
**160mm x 115mm**

**ASK FOR THE BRAND**

[shoprider.com.au](http://shoprider.com.au)