

# SHOPRIDER®

*mobility*

**10 YEAR**  
WARRANTY  
terms and conditions apply

**BLUE**  
888SEB

**SILVER**  
888SES

**RED**  
888SE

# EXPLORER

## NEW FEATURES

- Gas assist adjustable tiller
- USB port
- Comfortable high back contoured seat with swivel function
- Updated basket with carry handle
- Tie down points on rear of scooter
- Adjustable rear vision mirrors
- Rear adjustable suspension
- Height adjustable seat with depth adjustment away from tiller
- Width, height and angle adjustable flip up padded armrests
- Height adjustable headrest
- Angle adjustable backrest
- Padded loop handlebar
- Fingertip operated controls
- Forward / Reverse drive settings
- Variable speed settings
- Battery power indicator
- Horn
- LED Headlight
- Rear lights
- Rear wheel drive
- Pneumatic tyres
- Electromagnetic brake

NEW  
DASH



NEW  
MIRRORS



NEW  
BASKET

NEW  
STORAGE  
POCKET



**NEW**  
COLOUR OPTION!

Please contact **SHOPRIDER** for more information on each product

**Flip over for Specifications**

ASK FOR THE BRAND

[shoprider.com.au](http://shoprider.com.au)

**IMPORTANT NOTE:**  
It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

**SHOPRIDER** strives to provide the highest quality products possible that are subject to change without notice.

# SHOPRIDER®

*mobility*



NEW  
HEADLIGHT

NEW  
NUDGE BAR

NEW  
RIMS

# EXPLORER

## SPECIFICATIONS

- Overall: **1280mm (L) x 600mm (W) x 1070mm** (H at lowest settings)
- Seat: **460 - 530mm (W int.) x 445mm (D)**  
**650mm / 680mm / 695mm / 710mm / 725mm (H)**
- Backrest: **460mm (W) x 500mm (H) x 0 - 100°**  
(Angle range)
- Headrest: **255mm (W) x 180mm (H) x 630 to 750mm** (H from seat - 7 settings at 20mm increments)
- Armrests: **300mm (L) x 85mm (W) x 240 to 315mm** (H from seat - 4 settings at 25mm increments) **0 - 90°**  
(Angle range)
- Tiller: **460mm (W) x 900mm (H) x 30° forwards to 15° towards seat** (Angle range)
- Turning Circle: **1500mm**
- Speed: **9km/h**
- Maximum range: **Up to 32km**
- Safe climb angle: **9°**
- Maximum climb angle: **9°**
- Motor: **0.75Hp high torque motor (2 pole)**
- Batteries: **2 x 12V / 50 Amp (supplied)**
- Charger: **8 Amp (Supplied)**
- Tyres (Pneumatic): **Front / Rear: Ø260 x 85mm (10")**
- Product Weight: **99kg**
- Maximum User Weight: **136kg (on flat ground)**



0° - 136kg  
6° - 113kg  
9° - 91kg

P: **08 9248 4180**  
E: **orders@for-de.com.au**

ASK FOR THE BRAND

**shoprider.com.au**