

SHOPRIDER®

mobility

10 YEAR WARRANTY

terms and conditions apply

RED

STREAMER

888WA

FEATURES

- Easily maneuverable
- High back adjustable seat
- Sleek body design
- Adjustable arm rests and headrest
- Controller interchangeable both arms
- Free wheel device

Options

- Leg Rest Adaptors: **L: (770111A) R: (770112A)**
- Elevating leg rests: **L: (770113) R: (770114)**
- Swing away leg rests: **(770115)**
- Stump Support: **L: (SE-SSNL) R: (SE-SSNR)**

SPECIFICATIONS

- Length: **670mm** (front castors to rear tipper)
- Width: **580mm** (outside of rear wheels)
- Height: **1170mm** (ground to headrest)
- Seat Sizing: W: **490mm** x D: **480mm** x H: **560 - 610mm** (at highest setting)
- Adjustment range in armrests: W: **510 - 610mm**
- Battery Capacity (minimum amp): **2 x 36Amp**
- Maximum range: **up to 18km**
- Maximum load (single user): **136kg**
- Charger off board: **5Amp**
- Motor size: **0.8Hp x2**
- Maximum speed: **7.5km/h**
- Safe climbing angle: **9°**
- Max. climbing angle: **9°**
- Turning radius: **600mm**
- Wheel size: Front: **150mm** Back: **260mm**

Please contact **SHOPRIDER** for more information on each product

P: 08 9248 4180

E: orders@for-de.com.au

1 Business Way Malaga Western Australia 6090
Part of the **For-de Group Pty Ltd**

HEIGHT ADJUSTABLE
HEADREST

FOLDING
BACKREST

FLIP UP
ARM RESTS
(WIDTH /HEIGHT)

FLIP UP
FOOTREST

0° - 136kg
6° - 113kg
9° - 75kg

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

ASK FOR THE BRAND

shoprider.com.au