

**SHOPRIDER**<sup>®</sup>  
*Mobility*

**10 YEAR**  
**WARRANTY**  
terms and conditions apply

**RED**  
 889AER

**BLUE**  
 889AEB

**GREY**  
 889AEG

# ROCKY 8

## FEATURES

- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Headlight (LED)
- Rear lights
- LCD display
- Indicators / Hazard lights
- Large battery capacity
- Battery power indicator
- Horn
- USB power outlet
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket / storage pocket



0° - 225kg  
 6° - 150kg  
 9° - 125kg  
 12° - 100kg

## SPECIFICATIONS

- Overall: **1590mm (L) x 730mm (W) x 1285mm (H)**
- Turning Circle: **1860mm**
- Speed: **10km/h**
- Maximum range: **Up to 54km**
- Safe climb angle: **12°**
- Maximum climb angle: **12°**
- Motor: **2.0Hp high torque motor (4 Pole)**
- Batteries: **2 x 60 Amp minimum (2 x 75 amp supplied)**
- Charger: **10 Amp**
- Tyres (Pneumatic): **Ø330mm (80/80-8) - Front + Rear**
- Product Weight: **180kg - with batteries**
- Maximum User Weight: **225kg - on flat ground**

Please contact **SHOPRIDER**  
 for more information on this product

**P: 08 9248 4180**  
**E: orders@for-de.com.au**

**ASK FOR THE BRAND**

**shoprider.com.au**

**IMPORTANT NOTE:**  
 It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

**SHOPRIDER** strives to provide the highest quality products possible that are subject to change without notice.