

CHAMPAGNE 889RGC



X - GBC
Golf Bag Carrier
SOLD SEPERATELY

0° - 225kg
6° - 150kg
9° - 125kg
12° - 100kg

FEATURES

- Sliding / rotating / height adjustable seat
- High back contoured seat with headrest
- Lift up & adjustable width arm rests
- Steering tiller with 5 angle setting
- Straight handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + rear (adjustable) suspension
- Twin headlights (LED)
- Rear lights
- Indicators / Hazard lights
- Battery power indicator
- Horn
- Removable storage basket

SPECIFICATIONS

Overall: **1430mm (L) x 750mm (W) x 1300mm (H)**

Turning Circle: **1800mm**

Speed: **10km/h (up to 15km/h programmable)**

Maximum range: **Up to 45km**

Safe climb angle: **9°**

Maximum climb angle: **12°**

Motor: **2.0Hp high torque motor (4 Pole)**

Batteries: **2 x 55 amp minimum (2 x 75 amp supplied)**

Charger: **10 amp**

Tyres (Pneumatic): **Ø330 x 120mm (13 x 5.00-6)**

- Front + Rear

Product Weight: **166kg** - with batteries

Maximum User Weight: **225kg** - on flat ground

P: **08 9248 4180**

E: **sales@shoprider.com.au**

ASK FOR THE BRAND

shoprider.com.au