



**10 YEAR**  
WARRANTY  
terms and conditions apply

BLUE  
889WB

# ROCKY 12

**WIDE SEAT**  
OPTION ON ENQUIRY



## FEATURES

- High back contoured seat with headrest
- Wide seat option available
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Emergency hand brake
- Twin headlights (LED)
- Rear lights with brake-lights
- Full LCD display
- Automotive style Indicators / Hazard lights
- Large battery capacity 2 x 100Amps  
XT100 Premium included
- Battery power indicator
- USB power outlet
- Sliding seat
- Lift up & adjustable width arm rests
- Removable shopping basket
- 2 seat options

## SPECIFICATIONS

- Overall: 1740mm (L) x 750mm (W) x 1310mm (H)
- Std. Seat: 560mm (W) x 540mm (D) - (22" x 21")
- Wide Seat: 820mm (W) x 450mm (D) - (32" x 18")
- Turning Circle: 2705mm
- Speed: 10km/h
- Maximum range: Up to 50km
- Safe climb angle: 9°
- Maximum climb angle: 9°
- Motor: 2.0Hp high torque motor (4 Pole)
- Batteries: 2 x 100Ah (XT100G supplied)
- Charger: 8 Amp
- Tyres (Pneumatic): Ø410mm (100/80/10) - Front + Rear
- Product Weight: 186kg with batteries
- Maximum User Weight: 350kg - on flat ground

Please contact **SHOPRIDER** for more information on this product

P: **08 9248 4180**  
E: **orders@for-de.com.au**

**ASK FOR THE BRAND**

**shoprider.com.au**

**IMPORTANT NOTE:**  
It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

**SHOPRIDER** strives to provide the highest quality products possible that are subject to change without notice.