

SHOPPRIDER[®]
mobility

10 YEAR
WARRANTY
terms and conditions apply

BLUE
889XLB

RED
889XLR

YELLOW
889XLY

NEW
COLOUR

ROCKY 4

FEATURES

- NEW contoured high back swivel seat
- Sliding / rotating / height adjustable seat
- Angle adjustable backrest
- Height adjustable headrest
- Width / angle adjustable flip up armrests
- User friendly controls
- S drive control system
- Padded loop handlebar
- Fingertip operated forward / reverse drive
- Variable speed settings
- Battery power indicator
- A and C USB ports
- Angle adjustable tiller (gas assist)
- Tiller storage pocket
- Adjustable rear vision mirrors
- Removable basket with carry handle
- Horn
- LED twin headlights
- Indicators / hazard lights
- Tail / brake lights
- Rear wheel drive
- Front and rear adjustable suspension
- Pneumatic tyres
- Fluted rims
- Electromagnetic brake
- Front nudge bar



Flip over for specifications...

ASK FOR THE BRAND

shoprider.com.au

IMPORTANT NOTE:
It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

SHOPPRIDER strives to provide the highest quality products possible that are subject to change without notice.


SHOPRIDER®

mobility

ROCKY 4

SPECIFICATIONS

- Overall size (minimum settings): **1440mm (L) x 700mm (W) x 1300mm (H)**
- Turning Circle: **2850mm**
- Speed: **10km/h**
- Maximum range: **Up to 40km**
- Safe climb angle: **9°**
- Maximum climb angle: **12°**
- Motor: **2.0Hp high torque motor (4 pole)**
- Batteries: **2 x 75 Amp minimum (supplied)**
- Charger: **10 Amp**
- Front / rear tyres (Pneumatic): **Ø330 x 100mm / (4.00-5)**
- Product Weight: **146kg (with batteries)**
- Maximum User Weight: **225kg (on flat ground)**



Please contact **SHOPRIDER** for more information if required

P: 08 9248 4180

E: orders@for-de.com.au

ASK FOR THE BRAND

shoprider.com.au