



**10 YEAR**  
WARRANTY  
terms and conditions apply

BLUE  
889XLB

RED  
889XLR

# ROCKY 4

## FEATURES

- NEW contoured high back swivel seat
- Sliding / rotating / height adjustable seat
- Angle adjustable backrest
- Height adjustable headrest
- Width / angle adjustable flip up armrests
- User friendly controls
- S drive control system
- Padded loop handlebar
- Fingertip operated forward / reverse drive
- Variable speed settings
- Battery power indicator
- A and C USB ports
- Angle adjustable tiller (gas assist)
- Tiller storage pocket
- Adjustable rear vision mirrors
- Removable basket with carry handle
- Horn
- LED twin headlights
- Indicators / hazard lights
- Tail / brake lights
- Rear wheel drive
- Front and rear adjustable suspension
- Pneumatic tyres
- Fluted rims
- Electromagnetic brake
- Front nudge bar

**2024**  
MODEL UPGRADES



**Flip over** for specifications...

**ASK FOR THE BRAND**

**shoprider.com.au**

**IMPORTANT NOTE:**  
It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

SHOPRIDER strives to provide the highest quality products possible that are subject to change without notice.



# SHOPRIDER®

*mobility*

**2024**  
MODEL UPGRADES



NEW  
TILLER  
WITH  
STORAGE  
POCKET


NEW  
FLUTED RIMS

NEW  
FRONT END



NEW  
SEAT

NEW  
REAR END



INDICATORS


NEW  
USER  
FRIENDLY  
CONTROLS



# ROCKY 4

## SPECIFICATIONS

- Overall size (minimum settings):  
**1440mm (L) x 700mm (W) x 1300mm (H)**
- Turning Circle: **1800mm**
- Speed: **10km/h**
- Maximum range: **Up to 40km**
- Safe climb angle: **9°**
- Maximum climb angle: **12°**
- Motor: **2.0Hp high torque motor (4 pole)**
- Batteries: **2 x 75 Amp minimum (supplied)**
- Charger: **8 Amp**
- Front / rear tyres (Pneumatic):  
**Ø330 x 100mm / (4.00-5)**
- Product Weight: **146kg (with batteries)**
- Maximum User Weight:  
**225kg (on flat ground)**



0° - 225kg  
6° - 180kg  
9° - 159kg  
12° - 136kg

Please contact **SHOPRIDER**  
for more information if required

**P: 08 9248 4180**

**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**

*ASK FOR THE BRAND*

**[shoprider.com.au](http://shoprider.com.au)**