



888WNLM

BURGUNDY RED



FEATURES

- High back contoured seat with headrest
- Height / angle adjustable seat
- Lift up & adjustable width / height arm rests
- Fold up adjustable height / depth footrest
- Mid wheel drive
- Finger control operation
- High / Low speed settings
- Battery power indicator
- Horn
- Ideal for indoor and outdoor use on even ground
- Tie down points for transport
- Attachment options available:
 Swing away legrests
 Elevating swing away legrests
 Stump support

SPECIFICATIONS

Overall: 960mm (L) x 620mm (W) x 1160mm (H)

Turning Circle: 990mm Speed: 8km/h

Maximum range: Up to 26km

Safe climb angle: 8° Maximum climb angle: 12°

Motor: 2 x 0.85Hp high torque motors (2 Pole) Batteries: 2 x 36 Amp minimum (supplied)

Charger: 5 Amp

0° - 136kg 6° - 113kg

9° - 91kg

12° - 75kg

Tyres: Front: Ø180mm / Middle: Ø260mm / Rear: Ø180mm

Product Weight: 87kg - with batteries

Maximum User Weight: 136kg - on flat ground