



RED 9LSP

BLUE 9LSPB


GREY 9LSPG

ALLROUNDER

NEW FEATURES

- Easy access control
- User friendly tiller / dash
- USB charger
- Accessory / storage pocket
- Stylish fluted rims
- Double nudge bar
- LED indicators
- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear suspension
- LED Headlight
- Rear lights
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket





SPECIFICATIONS

- Overall: 1300mm (L) x 640mm (W) x 1200mm (H)
- Turning Circle: **1500mm**
- Speed: 10km/h
- Maximum range: **Up to 40km**
- Safe climb angle: 9°
- Maximum climb angle: 12°
- Motor: 1.5Hp high torque motor (4 pole)
- Batteries: 2 x 50 Amp minimum (supplied)
- Charger: 8 Amp
- Tyres (Pneumatic): Ø330 x 100mm (4.00-5) Front + Rear
- Product Weight: 99kg with batteries
- Maximum User Weight: 159kg on flat ground

Please contact **SHOPRIDER** for more information on each product

P: **08 9248 4180** E: **orders@for-de.com.au** ess Way Malaga Western Australia 6090

1 Business Way Malaga Western Australia 6090 Part of the **For-de Group Pty Ltd**

ASK FOR THE BRAND

shoprider.com.au