

Shoprider COMO

ASK FOR THE BRAND **SHOPRIDER**



SHOPRIDER[®]
www.shoprider.com.au
mobility

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER®
www.shoprider.com.au *mobility*

COMO

The SHOPRIDER™ **COMO** portable powerchair is a well featured powerchair that offers the user the ability to easily fold and transport by using quick disconnect levers in seconds. The **COMO** features a comfortable mid back seat, adjustable height armrests, easy to use joy stick control and when combined with the powerful motors it will deliver outstanding safety & comfort.

SPECIFICATIONS

Length	725mm
Width	545mm
Height	1180mm
Battery Capacity	12 Amp
#Maximum Range	up to 16km
Maximum Load Single User*	113Kgs
Charger on board	2Amp
Motor Size	0.50hp
Maximum Speed	5-6Kph
*Safe Climbing Angle	6 Degrees
*Max Climbing Angle	6 Degrees
Turning Radius	685mm
Wheel Size:	
- Front Caster	2 x 155mm
- Rear Caster	2 x 2 mm

FEATURES

12amp SLA Battery
Mid Back Captain Seat
Variable seat height
Adjustable arm rests
Controller interchangeable
Free wheel device
Easy Dismantle for Transport
Total Weight 47.4kg
Largest weight
18.4kg in batts front Chassis


Driver weight may exceed weight of unit.
Reduce speed when turning.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

COMO



SHOPRIDER®
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au