



10 YEAR
WARRANTY
terms and conditions apply

RED
G424M

COUGAR 10

FEATURES

- High back contoured seat with headrest
- Finger control operation
- High / Low speed settings
- Mid wheel drive
- Battery power indicator / Horn
- Height / angle adjustable seat
- Lift up & adjustable width / height arm rests
- Fold up adjustable height / depth footrest
- Ideal for indoor and outdoor use on even ground
- Tie down points for transport

Options

- Adaptor brackets: **77011A (L) x 77012A (R)** (required)
- Elevating swing away leg rests: **770113 (L) x 770114 (R)**
- Swing away leg rests: **770115**

SPECIFICATIONS

- Overall: **920mm (L) x 620mm (W) x 1120mm (H)**
- Turning Circle: **950mm**
- Speed: **7-8km/h**
- Maximum range: **Up to 40km**
- Safe climb angle: **8°**
- Maximum climb angle: **9°**
- Motor: **2 x 0.8Hp high torque motors (2 Pole)**
- Batteries: **2 x 50 Amp minimum (supplied)**
- Charger: **5 Amp**
- Tyres (Pneumatic): **Ø140 x 50mm (Front)**
Ø260 x 85mm (Middle)
Ø127 x 37mm (Rear)
- Product Weight: **65kg - with batteries**
- Maximum User Weight: **136kg - on flat ground**



0° - 136kg
6° - 113kg
9° - 91kg

Please contact **SHOPRIDER** for more information on this product

ASK FOR THE BRAND

P: **08 9248 4180**
E: **orders@for-de.com.au**

shoprider.com.au

IMPORTANT NOTE:
It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

SHOPRIDER strives to provide the highest quality products possible that are subject to change without notice.