

778HD BURGUNDY RED



FEATURES

- High back contoured seat with headrest
- 5 point adjustable front tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear suspension
- Headlights
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

SPECIFICATIONS

Overall: 1230mm (L) x 640mm (W) x 1250mm (H)
 Turning Circle: 1250mm
 Speed: 10km/h
 Maximum range: Up to 32km
 Safe climb angle: 9°
 Maximum climb angle: 9°
 Motor: 2.0Hp high torque motor (2 Pole)
 Batteries: 2 x 50 Amp minimum (supplied)
 Charger: 5 Amp
 Tyres (Flat free): Ø260mm Front / Ø330mm Rear
 Product Weight: 102kg - with batteries
 Maximum User Weight: 182kg - on flat ground

