


889D
BURGUNDY RED



0° - 225kg
6° - 180kg
9° - 159kg
12° - 136kg

FEATURES

- Comfortable large wide bench seat - 800mm
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Headlight
- Rear lights
- Indicators
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding seat
- Lift up & adjustable width arm rests
- Removable shopping basket

SPECIFICATIONS

Overall: 1480mm (L) x 910mm (W) x 980mm (H)
Turning Circle: 2550mm
Speed: 10km/h
Maximum range: Up to 40km
Safe climb angle: 9°
Maximum climb angle: 12°
Motor: 2.0Hp high torque motor (2 Pole)
Batteries: 2 x 75 Amp minimum (supplied)
Charger: 8 Amp
Tyres (Pneumatic): Ø330mm - Front + Rear
Product Weight: 159kg - with batteries
Maximum User Weight: 225kg - on flat ground

SINGLE PERSON USE ONLY