


**GK10**  
BURGUNDY RED

**GK10B**  
SHINING BLUE



0° - 136kg  
6° - 113kg  
9° - 91kg

## FEATURES

- Scooter pulls apart for storage and transport
- Comfortable mid back contoured seat
- Adjustable steering tiller
- Loop handle with finger operated controls
- Rear wheel drive
- Electromagnetic brake
- Battery power indicator
- Split battery box
- Headlights
- Removable rotating / height adjustable seat
- Lift up & adjustable width arm rests
- Lifting handles makes picking up the scooter easy
- Removable shopping basket

## SPECIFICATIONS

Overall: 1190mm (L) x 585mm (W) x 960mm (H)  
Turning Circle: 1330mm  
Speed: 9km/h  
Maximum range: Up to 27km  
Safe climb angle: 9°  
Maximum climb angle: 9°  
Motor: 1.2Hp high torque motor  
Batteries: 2 x 36 Amp minimum (supplied)  
Charger: 5 Amp  
Tyres (Pneumatic): Ø260mm - Front + Rear  
Maximum User Weight: 136kg - on flat ground

Product Weight: 78kg - with batteries  
Battery Weight: 11kg each  
Seat Weight: 16kg  
Rear Transaxle Weight: 17kg  
Main Chassis Weight: 23kg