



## FEATURES

GK9-3

- Scooter pulls apart for storage and transport
- Comfortable mid back contoured seat
- Fold down adjustable steering tiller
- Straight handle with finger operated controls
- Rear wheel drive
- Electromagnetic brake
- Battery power indicator
- Split battery box
- Horn
- Removable rotating / height adjustable seat
- Lift up & adjustable width arm rests
- Lifting handles makes picking up the scooter easy
- Removable shopping basket

## SPECIFICATIONS

Overall: 990mm (L) x 560mm (W) x 880mm (H)  
 Turning Circle: 950mm  
 Speed: 6km/h  
 Maximum range: Up to 20km  
 Safe climb angle: 9°  
 Maximum climb angle: 9°  
 Motor: 0.8Hp high torque motor (2 Pole)  
 Batteries: 2 x 21 Amp minimum (supplied)  
 Charger: 2 Amp  
 Tyres (Flat free): Ø230mm - Front + Rear  
 Maximum User Weight: 136kg - on flat ground

Product Weight: 50kg - with batteries  
 Battery Weight: 7kg each  
 Seat Weight: 8.3kg  
 Rear Transaxle Weight: 11.2kg  
 Main Chassis Weight: 16.5kg