



LITTLE RIPPER

PORTABLE MOBILITY SCOOTER

GK9 - RED **GK9** - SHINING BLUE

FEATURES

- Pull apart design for transport
- Mid back Fold down seat
- Rotating / Height adjustable seat
- Adjustable T Bar steering tiller
- Rear wheel drive
- Electromagnetic brake
- Battery power indicator
- Lift up adjustable width arm rests
- Removable shopping basket
- Front bumper bar

SPECIFICATIONS (All dimensions shown are nominal)

Overall: 1050mm (L) x 560mm (W) x 880mm (H)

Turning Circle: 1150mm

Speed: 6kmph

Maximum range: Up to 20km

Safe climb angle: 9° Maximum climb angle: 9° Motor: 0.8Hp motor

Batteries: 2 x 21 amp minimum

Charger: 4 amp

Tyres: Ø225mm flat free

Ground Clearance: 45mm (30mm - tipper wheel)

Product Weight: 56kg - with batteries

Maximum User Weight: 136kg - on flat ground



Scooter pulls apart into manageable pieces for transport:

- Main Chassis: 16.4kg - Rear Transaxle: 14.3kg

- Chair: 10.3kg

- Battery: 7.5kg each

Please contact **SHOPRIDER** if you require further information on this product.

P: 08 9248 4180

F: 08 9248 9452

E: orders@shoprider.com.au

SHOPRIDER - 1 Business way, Malaga, Western Australia 6090 Part of the **FOR-DE Group** Pty Ltd.

shoprider.com.au

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

Easy to pull apart and put together

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.