



889XLSBNRR
INJECTION RED

889XLSBNRRB
INJECTION BLUE



0° - 160kg
6° - 110kg
9° - 90kg
12° - 70kg

RAINRIDER

FEATURES

- Fully enclosed cabin with double sided doors
- Windscreen wipers
- Double pillow comfort high back seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Twin Headlights (LED)
- Rear lights
- Indicators / Hazard lights
- Large battery capacity
- Battery power indicator
- USB power outlet
- Sliding / Rotating / Height adjustable seat

SPECIFICATIONS

Overall: 1600mm (L) x 790mm (W) x 1620mm (H)
Turning Circle: 1800mm
Speed: 10km/h
Maximum range: Up to 40km
Safe climb angle: 9°
Maximum climb angle: 12°
Motor: 2.0Hp high torque motor
Batteries: 2 x 75 Amp minimum (supplied)
Charger: 8 Amp
Tyres (Pneumatic): Ø330mm - Front + Rear
Product Weight: 195kg - with batteries
Maximum User Weight: 160kg - on flat ground