

889RGB
 DARK BLUE




Optional
GOLF BAG CARRIER (SE-GBM)
 can be fitted to scooter
 (Not included - Sold separately)

FEATURES

- High back contoured seat with headrest
- Adjustable steering tiller
- Straight handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Twin Headlights (LED)
- Rear lights
- Indicators
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable storage basket

SPECIFICATIONS

Overall: 1430mm (L) x 750mm (W) x 1300mm (H)
 Turning Circle: 1800mm
 Speed: 10km/h
 Maximum range: Up to 45km
 Safe climb angle: 9°
 Maximum climb angle: 12°
 Motor: 2.0Hp high torque motor (4 Pole)
 Batteries: 2 x 55 Amp minimum (not supplied)
 Charger: 10 Amp
 Tyres (Pneumatic): Ø330 x 120mm - Front + Rear
 Product Weight: 146kg - with batteries
 Maximum User Weight: 225kg - on flat ground



0° - 225kg
 6° - 150kg
 9° - 125kg
 12° - 100kg