



NEW PRODUCT RELEASE!



LIFT
position

VIBRATION MASSAGE
handset



RECLINE
position



TOFFEE
(LC0201)



GREY
(LC0202)



BLACK
(LC0203)

VITTORIA

DUAL MOTOR LIFT CHAIR

LC0201 - TOFFEE | LC0202 - GREY | LC0203 - BLACK

FEATURES

- Lift / Rest / Recline / Trendelenburg armchair
- Allows the user to be easily raised to a standing position and offers the comfort of a full recline
- Dual motors for the ultimate in adjustability
- Inbuilt vibration massage function
- Added comfort of a pocket spring seat base
- Luxurious leather coverings
- Available in modern colours: Toffee / Grey / Black
- Quality German actuation with battery backup
- Warranty:
 - Frame: 3 years / Electrics: 2 years / Fabric: 2 years

SPECIFICATIONS

All dimensions shown are nominal

Overall: **905mm (W) x 970mm (D) x 1075mm (H)**
 Seat: **470mm (W) x 520mm (D) x 500mm (H)**
 Armrest: **160mm (H) x 500mm (L)**
 Backrest: **770mm (W) x 670mm (H): 93 - 185°**
 Legrest: **630mm (W) x 540mm (L): 90 - 185°**
 Distance from wall: **590mm**
 Product Weight: **55kg**
 Maximum User Weight: **160kg**

Please contact **REDGUM Brand** if you require further information

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090

redgumbrand.com.au

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.